




JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING

EVT24: Medical Device for Gastrointestinal Leaks



Alex Ren | Damian Cano | Jay Heymann | Thomas Miller

Problem Scenario


About 1 in 5 gastrointestinal leaks lead to death




Problem Description



Endoscopic Vacuum Therapy






Extra-luminal



Intra-luminal

Cavity




Vacuum Device


Healed Cavity

Current State

Vacuum Stents




Sponges



Takeaway: There is no FDA approved EVT device, leaving physicians to construct “makeshift” devices.

EVT Procedure



Design Requirements

Procedure Time
under an Hour



Remain in Body
for a Week




Mass
Production



Endoscopic
Deployment




Active Vacuum




FDA Approval




Device and Subsystems




Video + Animation




Improvements Over Endo/Eso Sponge




Use of
Overtube




Passage
From Mouth
to Nose



Post Endoscopic
Evaluation




Streamlined Deployment




Simultaneous
Endoscopic Evaluation

✓ Decreased Procedure Time



✓ Reduced Required Equipment



✓ Increased Patient Comfort



Gelatin Casing

0.5 mm wall thickness

Ø 6mm

50mm



**Simulated Pill Capsule
Formula**



Porcine Gelatin



Water Soluble

Gel Casing Manufacturing



Mix / Heat Gelatin Solution



Lower Mold Into Solution



Allow Casing to Cure



Remove Casing From Mold



Alternative Molds With Cone

Gelatin Casing Testing


Partial Dissolution



Complete Dissolution





Procedure Steps




Average Partial Dissolution Time [min:sec]	Average Complete Dissolution Time [min:sec]	Estimated Procedure Time
6:50	8:53	10:00

Sponge Design


- Open-cell, cylindrical, 45 pores per inch, polyurethane sponge
- Biggest risk is tissue ingrowth, which is mitigated with two different coatings
 - Polyethylene glycol (**PEG**)
 - Poly lactic co-glycolic acid (**PLGA**)




PEG-Coated Sponge



PEG Solution



No PEG



15% PEG

Images taken using scanning electron microscopy at 10x zoom.

Results: Inconclusive as there is no significant bacterial growth in either case.

Bioabsorbable Layer




Polyurethane Sponge Covered
in Bioabsorbable Layer



Bioabsorbable Layer
(PLGA Scaffold)

Dissolvability Testing



Results: Achieves 7 day minimum but rate of dissolution is slow.

Steps Forward




Requirements Closeout

Procedure Time
under an Hour



Remain in Body
for a Week



Mass
Production



Endoscopic
Deployment



Active Vacuum




FDA Approval




Budget and Unit Cost

PLGA-PU Hybrid Sponge
Unit Cost: **\$140.47**

PEG-Coated Sponge Unit
Cost: **\$8.23**



Conclusion





JOHNS HOPKINS

WHITING SCHOOL *of* ENGINEERING

Alex Ren | Damian Cano | Jay Heymann | Thomas Miller